Professional cleaning and disinfection technology



Practical and versatile: MEIKO racks and holders



Focus on users: MEIKO cleaning and disinfection appliances





MEIKO cleaning and disinfection appliances not only meet the highest standards of safety, cleanliness and hygiene – they are also designed to be highly efficient and easy to use.

Our cleaning and disinfection appliances offer an impressively wide range of functions and features. One example is the MEIKO range of wash racks and holders which make it easier than ever to get the results you want.

An all-purpose, universal rack is included as standard with all MEIKO cleaning and disinfection appliances. This rack is designed to accommodate all the care utensils typically encountered in a nursing environment, including bedpans of various shapes and sizes, lids, commode pots, kidney bowls, and numerous other items.

MEIKO can also provide custom racks and holders tailored to your specific requirements. All our racks are designed to be functional and practical – and switching from one rack to another couldn't be easier!

Our flexible quick-change system works with all MEIKO racks and holders, true to our policy of offering the ultimate in user-friendliness.



Focus on benefits: MEIKO universal rack

Exceptional benefits that make life easier for care personnel

- Universal rack available in two different versions (SUH 1 or SUH 2)
- Simply place the utensils in the rack – no need to empty them first
- Open door also acts as a drip tray
- As soon as you close the door the utensils are emptied in a sealed environment
- Rack designed to accommodate all the standard care utensils used on wards
- Switching from one rack to another is quick and easy, offering greater flexibility in everyday work situations
- Optional accessory: Basket inserts for small items which fit neatly into the universal rack

Exceptional benefits that promote excellence in nursing care

- Easier handling of care utensils
- Superb standards of hygiene
- Reliable, user-friendly operation

MEIKO technology for clean solutions.



Universal rack (SUH 1)



Kidney bowl / urine bottle (SUH 1 or SUH 2)



Urine jug / urine bottle (SUH 1 or SUH 2)





Wash basin (SUH 1 or SUH 2)



3 urine bottles (SUH 1 or SUH 2)



Fine-mesh basket insert to keep small items separate (optional accessory) (SUH1 or SUH2)



Bedpan with lid / urine bottle (SUH 1 or SUH 2)



Commode pot with lid / urine bottle (SUH 1 or SUH 2)



Child's bedpan with handle and lid (SUH 1)



Oval bedpan / urine bottle (SUH 2)



Slipper-style bedpan with lid / urine bottle (SUH 2)

Focus on your needs: MEIKO custom racks and holders

We specialise in tailor-made solutions

Whatever care utensils you use, MEIKO can provide you with custom racks and holders to empty and rinse them in a sealed environment. Our flexible quick-change system works with all MEIKO racks and holders, offering you the ultimate in user-friendliness.

The right adapter for every application

For special utensils that are only used every now and then, MEIKO offers custom adapters that can be inserted without having to remove the universal rack.



Quick-change rack for child's potty



Quick-change rack for Redon bottles



Quick-change rack for 4 urine bottles



Quick-change holder for suction bottle and lid



Rack with adapter for suction bottles



Quick-change holder for collection utensils and lids



Rack with adapter for urine jugs



Care utensil wash racks for worldwide use

MEIKO cleaning and disinfection appliances are used all over the world and are adapted to meet the cleaning requirements of each individual country. MEIKO provides custom racks and holders for a wide range of countries designed to accommodate the care utensils used in that particular region of the world.

Six examples of our country-specific racks are shown here.

















MEIKO Maschinenbau GmbH & Co. KG Englerstr. 3 77652 Offenburg Germany Tel.: +49(0)781 203-0 Fax: +49(0)781 203-1229 E-Mail: info@meiko.de www.meiko.de

740.061.01.04.16/GB/KD/web

We reserve the right to amend specifications as part of our product improvement process.